



ANKLE & FOOT

POST OPERATIVE SHOE	29	LOWER LEG CORK SPLINT	38
POST OPERATIVE HIGH TOP	30	PRE-CUT MOULDED SPLINT	38
ALL PURPOSE BOOT	31	FIXED WALKER BOOT	39
MED-SURGE SHOE	32	ROM WALKER BOOT	40
ORTHO WEDGE SHOE	33	MEMORY ANKLE BRACE	41
WEDGE HEEL SHOE	34	ECO WALKER	41
PEG ASSIST SYSTEM	35	AIR/GEL ANKLE BRACE	42
PEG ASSIST INSOLE	35	EXO ANKLE BRACE	42
SWEDISH AFO	36	SILICONE HEEL CUP	43
PLANTAR FACIITIS SUPPORT	36	DARCO TOE ALIGNMENT SPLINT	43
NIGHT SPLINT	37	ACHILLES HEEL PADS	43
FX NIGHT BOOT	37		

BENEFOOT™ 'ORIGINAL' POST-OP MEDICAL SHOE

The **BeneFoot™ 'Original'** can be used for post-op procedures or any indications that prevent the use of normal shoes and has been created using black nylon mesh to provide shoe like features for continued comfort and support throughout the day.



FEATURES AND BENEFITS:

- Made with a secure polypro shank built into the outsole creating a fully rigid 'shoe-like' sole. Allows for better flexion control.
- Round toe configuration provides protection and comfort.
- Closure system enables secure fit with added dressings.
- Dual density insole - Comfortable, allows for insole modification.
- Flared TPR outsole with rocker bottom, aids in ambulation and providing increased stability.
- Padded collar & top line to reduce heel slippage.
- Patient wellbeing - Better aesthetics, looks more like a shoe.
- Reinforced heel counter - For better medial/lateral stability.
- Quality construction - 100% Satisfaction guarantee



Aid in abulation and provides better stability



Padded collar & topline



Reinforced heelcounter

INDICATIONS FOR USE:

- For recovery after post-operative and tissue procedures.
- For use with any size leg casts.
- To accommodate large dressings or swellings.
- Diabetic foot indications.
- Any situation where size of foot prohibits normal footwear.

Post-Op Medical Shoe Sizing Information:

Code	Size	UK Female	UK Male	EU
POPD	Paediatric			
POWS	Women's Small	2.5 - 4.5		35 - 37.5
POWM	Women's Medium	5 - 6.5		38 - 39.5
POWL	Women's Large	7 - 8.5		40 - 43
POMS	Men's Small		5 - 7	37.5 - 40
POMM	Men's Medium		7.5 - 9	41 - 43
POML	Men's Large		9.5 - 11	44 - 46
POMXL	Men's X-Large		11.5 - 13	46.5 - 48.5

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

BENEFOOT™ 'HIGH TOP' MEDICAL SHOE

The BeneFoot™ 'High Top' has all the benefits of the Post-Op Original, plus added ankle support provided by the distinctive high top design creating less slippage and greater stability.

FEATURES AND BENEFITS:

- Eliminates heel slippage and forward foot movement.
- Round toe configuration provides protection and comfort.
- Closure system enables secure fit with added dressings.
- Made with a secure polypro shank creating a fully rigid 'shoe-like' sole.
- Flared Black TPR Outsole - Provides better stability.
- Reinforced Heel Counter - For better Medial/Lateral stability.
- Patient wellbeing - Better aesthetics, looks more like a shoe.
- Quality construction - 100% Satisfaction guarantee
- Eliminates Pressure / Compression on the toes due to Nylon Mesh Upper ends behind MTP joint.

INDICATIONS FOR USE:

- Plantar Wounds.
- After post-operative and tissue procedures.
- For use with any size leg casts.
- To accommodate large dressings or swellings.
- Diabetic foot indications.
- Any situation where size of foot prohibits normal footwear.

High Top Medical Shoe Sizing Information:

Code	Size	UK Female	UK Male	EU
POHTWS	Women's Small	2.5 - 4.5		35 - 37.5
POHTWM	Women's Medium	5 - 6.5		38 - 39.5
POHTWL	Women's Large	7 - 8.5		40 - 43
POHTMS	Men's Small		5 - 7	37.5 - 40
POHTMM	Men's Medium		7.5 - 9	41 - 43
POHTML	Men's Large		9.5 - 11	44 - 46
POHTMXL	Men's X-Large		11.5 - 13	46.5 - 48.5

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE



Better aesthetics to make it look more like a shoe



Better aesthetics to make it look more like a shoe



DARCO ALL PURPOSE BOOT

The **Darco All Purpose Boot** is a boot for all seasons and its closed toe protection allows dressings to be kept clean and toes dry, providing more protection for post-op procedures and leg casts.

FEATURES AND BENEFITS:

- **Square toe design acts as a bumper when K-wires are present and provides additional room and comfort for the patient.**
- **Wider sole and deep opening gives more room to accommodate large dressings, bandages and casts.**
- **Wide ankle strap eliminates heel slippage as well as secures the foot or cast in the shoe.**
- **Rocker bottom encourages smooth roll through gait pattern and reduces plantar pressure on the heel and forefoot**
- **Strapless forefoot closure provides uniform security without the pressure of buckles or straps.**
- **Insole can be removed for improved accommodation and in conjunction with peg-assist system.**
- **Quality construction - 100% Satisfaction guarantee**

INDICATIONS FOR USE:

- **Post-operative and tissue procedures.**
- **For use with any size leg casts.**
- **To accommodate large dressings or swellings.**
- **Diabetic foot indications.**

All Purpose Boot Sizing Information:

Code	Size	UK Women	UK Male	EU
APB-1	X-Small	2-5		34-37.5
APB-2	Small	5.5-8	5.5-7.5	38-41.5
APB-3	Medium	8.5-11	8-9.5	42-44
APB-4	Large	11+	10-11.5	44.5-46.5
APB-5	X-LLarge		12-13.5	47-49

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE



Square toe design



Wide ankle strap





DARCO MED SURG SHOE

The **Darco Med-Surg Shoe's** is a fully adjustable and versatile. Its strapless forefoot closure eliminates buckle pressure on the forefoot and can expand to accommodate dressings up to 17" in circumference.

FEATURES AND BENEFITS:

- Square toe design acts as a bumper to provide additional protection when K-wires are present.
- Rocker sole reduces plantar pressure on the forefoot and heel.
- New ankle strap eliminates heel slippage
- Zoned outsole provides traction under the midfoot
- MetaShank™ Protection provides more rigid control under the metatarsal heads
- Removable insole fully customizable with over twice the padding
- Quality construction - 100% Satisfaction guarantee

INDICATIONS FOR USE:

- Post-operative and tissue procedures.
- To accommodate large dressings or swellings.
- Diabetic foot indications.
- As a platform for the peg-assist insole.

Med-Surg Shoe Sizing Information:

Code	Size	UK Women	UK Male	EU
Med-surg 1	Paediatric			
Med-surg 2	Women's Small	2-5		
Med-surg 3	Women's Medium	5.5-8		
Med-surg 4	Women's Large	8.5-11		34-37.5
Med-surg 5	Men's Small	11+	5.5-7.5	38-41.5
Med-surg 6	Men's Medium		8-9.5	42-44
Med-surg 7	Men's Large		10-11.5	44.5-46.5
Med-surg 8	Men's X-Large		12-13.5	47-49

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE



BENEFOOT™ WEDGE SHOE

The **Benefoot™ Wedge** has been clinically proven to effectively reduce weight-bearing pressure from the forefoot to promote enhanced healing after surgery. The **Benefoot™ Ortho Wedge** shoe provides a lower degree of wedge for enhanced patient comfort whilst continuously reducing forefoot pressure.

Ortho Wedge Heel Shoe Sizing Information:

Code	UK Women	UK Male	EU
WSXS	2-5		34-37.5
WSS	5.5-8	5.5-7.5	38-41.5
WSM	8.5-11	8-9.5	42-44
WSL	11+	10-11.5	44.5-46.5
WSXL		12-13.5	47-49

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

FEATURES AND BENEFITS:

- Square toe design acts as a bumper to provide additional protection when K-wires are present.
- Advanced tread system increases traction where needed most.
- 15 Degree Dorsiflection to the midfoot and heel.
- Wide base for extra stability.
- Removable forefoot closure provides security whilst eliminating buckle pressure.
- Insole can be removed for improved accommodation and in conjunction with peg-assist system.
- Quality construction - 100% Satisfaction guarantee

INDICATIONS FOR USE:

- Post-surgical healing.
- Forefoot trauma.
- Wound or ulcerations present on the base of the foot.
- Any situations requiring a pressure offload of the forefoot.
- Diabetic foot indications.



DARCO WEDGE HEEL SHOE

The **Darco Wedge Heel Shoe** has been designed to offload the heel by shifting weight to the midfoot to promote faster healing after surgery, trauma, or when wounds are present on the heel.

FEATURES AND BENEFITS:

- **Reduced heel pressure.**
- **Square toe design acts as bumper to protect toes.**
- **Modified rocker bottom encourages a smooth toe-off while maintaining a stable mid-foot base.**
- **Unique forefoot closure provides uniform security without the pressure points of buckles and strap.**
- **Ankle strap seats the foot firmly in the sole and eliminates heel friction.**
- **Cloth toe protector keeps toes warm and dressings clean and dry.**
- **Insole can be removed for improved accommodation and in conjunction with peg-assist system.**
- **Quality construction - 100% Satisfaction guarantee**

INDICATIONS FOR USE:

- **Rear Foot trauma.**
- **Wounds and ulcerations on the heel area.**
- **Post-surgical healing for either soft tissue or bone structure of the heel.**
- **Diabetic foot indications.**

Wedge Heel Shoe Sizing Information:

Code	UK Women	UK Male	EU
WHSXS	2-5		34-37.5
WHSS	5.5-8	5.5-7.5	38-41.5
WHSM	8.5-11	8-9.5	42-44
WHSL	11+	10-11.5	44.5-46.5
WHSXL		12-13.5	47-49

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

DARCO PEG-ASSIST SYSTEM + INSOLE

The **Darco Peg-Assist System** is built around a multi-purpose removable peg insole that effectively off-loads the plantar aspect of the foot. This system comes with a dedicated toe covering to protect the foot and bandages from dirt and cold. The Peg Assist Insole is also available separately.



FEATURES AND BENEFITS:

- Easily removable pegs that off-load pressure from any area of the plantar surface of the foot.
- Constructed of an 18mm combination of Plastazote®, Poron®, and EVA to provide the perfect combination of shock absorption, off-loading and mould ability.
- Pressure is reduced allowing enhanced patient mobility.
- The PAS comes with a self-adhesive Peg Support Board that keeps pegs from collapsing and the off-loaded area intact.

INDICATIONS FOR USE:

- Available as an insole for most BeneFoot orthopaedic footwear including All Purpose Boot, Medsurg, Ortho Wedge and Heel Wedge.



Peg Assist System Sizing Information:

Code	Size	UK Women	UK Male	EU
PAS-1	X-Small	2-5		34-37.5
PAS-2	Small	5.5-8	5.5-7.5	38-41.5
PAS-3	Medium	8.5-11	8-9.5	42-44
PAS-4	Large	11+	10-11.5	44.5-46.5
PAS-5	X-Large		12-13.5	47-49

Peg Assist Insole Sizing Information:

Code	Size	UK Women	UK Male	EU
PASTC-1	Paediatric			
PASTC-2	Women's Small	2-5		
PASTC-3	Women's Medium	5.5-8		
PASTC-4	Women's Large	8.5-11		34-37.5
PASTC-5	Men's Small	11+	5.5-7.5	38-41.5
PASTC-6	Men's Medium		8-9.5	42-44
PASTC-7	Men's Large		10-11.5	44.5-46.5
PASTC-8	Men's X-Large		12-13.5	47-49

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

BENECARE™ SWEDISH AFO

The **BeneCare™ Swedish AFO** is a comfortable, low profile prefabricated leaf spring design to correct flaccid foot drop, especially on post-stroke patients. Injection moulded polypropylene allows for different thicknesses throughout the orthosis, thicker on the vertical aspect for rigidity, thinner footplate for ease of trimming, comfort and flexibility.

FEATURES AND BENEFITS:

- Injection-moulded for a high quality finish
- Easily trimmed footplate
- May be heat modified for minor fitting adjustments

INDICATIONS FOR USE:

- Drop foot
- Mild drop foot secondary to other neurological pathologies

Swedish AFO Sizing Information:

Left	Right	Size	UK	EU
SAFOLS	SAFORS	Small	1-3.5	33-36
SAFOLM	SAFORM	Medium	4-7	37-40
SAFOLL	SAFORL	Large	7.5-9.5	41-44
SAFOLXL	SAFORXL	Xlarge	10-14	45-48

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE



36



BENECARE™ PFS SUPPORT

The **BeneCare™ Plantar Fasciitis Support** provides continuous elastic tension and pressure to reduce the pain associated with inflammation due to plantar fasciitis. The PFS reduces pain, allowing the patient to strengthen the muscles of the foot, arch and lower leg, as well as stretch the Achilles tendon to provide relief and prevent future discomfort from plantar fasciitis.

FEATURES AND BENEFITS:

- Simple and inexpensive way to reduce inflammation.
- Elastic tension helps relieve pain and promote muscle development.
- Anatomically designed to provide comfort and control. Available in two adjustable sizes.
- Maintains tension continuously during wear, allowing the patient to exercise.
- Fast, easy 15-second application.
- Breathable stiff foam, wickable fabric draws perspiration away from skin.

PFS Support Sizing Information:

Code	Description	UK	EU
FPFS-SL	Small Left	3 to 5	35 - 38
FPFS-SR	Small Right	3 to 5	35 - 38
FPFS-LL	Large Left	6 to 14	38.5 - 49
FPFS-LR	Large Right	6 to 14	38.5 - 49

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

INDICATIONS FOR USE:

- Mild to moderate plantar fasciitis inflammations.

BENECARE™ NIGHT SPLINT

The **BeneCare™ Night Splint** is a lightweight and comfortable alternative to serial casting, for the treatment of plantar fasciitis, Achilles tendonitis and other lower extremity overuse injuries. Constructed from a strong durable plastic, the **BeneCare™ Night Splint** encases the foot and ankle in a Lycra® lined padded rigid shell, which helps to ensure total immobilization while maintaining comfort across the extensor tendons.

FEATURES AND BENEFITS:

- Adjustable flexion straps provide gentle stretching of the plantar fascia and Achilles tendon.
- Soft, Lycra® lined foam interior padding for optimal comfort.
- Breathable to eliminate heat build-up.
- Washable to aid with personal hygiene.
- Fully adjustable straps to maintain the optimum position and patient comfort.
- User-friendly, centre-release buckles simplify application.
- Slip-Resistant Tread.

Night Splint Sizing Information:

Code	Size	UK Male	UK Female	EU Size
PFSS	Small	Up to 6	Up to 5	Up to 39
PFSM	Medium	7-9.5	5.5-8	38-43
PFSL	Large	10+	8.5+	43.5+

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

INDICATIONS FOR USE:

- Plantar fasciitis
- Achilles tendonitis repair
- Lower extremity overuse injuries
- Prevention of drop foot



37

DARCO™ FX NIGHT SPLINT

The **Darco FX Night Boot** is designed to alleviate heel pain associated with plantar fasciitis and Achilles tendonitis. Two 1/2 inch EVA insoles are available to further the dorsiflex position of the foot.

FEATURES + BENEFITS:

- Non-slip treads on the sole.
- The wide planter aspect of the splint will provide excellent support for the indication.

INDICATIONS FOR USE:

- Plantar Fasciitis
- Heel Pain
- Tendinitis in heel area
- Prophylaxis
- For immobilisation

FX Night Splint Sizing Information:

Code	Size	UK Size	EU Size
FXNBS	Small	1-3.5	33-36
FXNBM	Medium	4-7	37-40
FXNBL	Large	7.5-9.5	41-44
FXNBXL	X-Large	10-14	45-48

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE



Best for

Tendon shortening, heel tendinitis and achilles tendon injuries

BENECARE™ LOWER LEG CORK SPLINT

The **UFS system** is an innovative, universal, ready to use splinting system made out of cork which is also the better substitute for Plaster of Paris (P.O.P.). As the system is an easy to use, fast and clean product with a high level of patient comfort it is a very versatile splinting system for both fracture clinics and A&E departments where working quick and efficiently with minimal mess will save time and money and improve patient care.

INSTRUCTIONS FOR USE:

- Select the appropriate size of splint. If needed, the splint can easily be adjusted with scissors to the desired measurements
- Slip stockinette over the effected body part and wrap additional padding over the stockinette. Padding should be used for protection of pressure points and bony prominences
- Fixate the splint using hook and loop straps
- The UFS-splint will automatically follow the contours of the extremities

Lower Leg Cork Splints Sizing Information:

Code	Description	Size
UFCS-03S	Low Leg	Adult Small
UFCS-03M	Low Leg	Adult Medium
UFCS-03L	Low Leg	Adult Large
UFCS-13S	Low Leg	Paediatric Small
UFCS-13M	Low Leg	Paediatric Medium
UFCS-13L	Low Leg	Paediatric Large

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE



BENECARE™ PRE-CUT MOULDED SPLINT

BeneCast™ Pre-Cut Moulded Leg Splints are an excellent way of applying strong durable splints. Simply wet the splint, by spraying or running under the tap and roll in a towel to absorb excess moisture. Apply directly to the limb using a Benecare cotton crepe bandage. Will mould perfectly to the contours of the limb and will reach functional strength in 30 mins.

FEATURES AND BENEFITS:

- Can be used quickly in emergency situations by ambulance services.
- 100% penetration of X-ray.
- 5-times lighter and stronger than Plaster Casting.
- Fast and efficient method to make a removable splint.
- Can be secured by either bandage or Velcro® tape.
- Pre-cut to its individual shape the splint requires minimal trimming or preparation.
- Excellent conformability for the perfect fit to the contours of the body.

Pre-Cut Moulded Leg Splints Sizing Information:

Code	Description	Size
PCMS010	Short Leg	Small
PCMS007	Short Leg	Large
PCMS006	Long Leg	Small

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE



BENECARE™ FIXED WALKER

A lightweight **Fixed Walker** providing support and comfort to patients.

FEATURES AND BENEFITS:

- Hook & Loop closure provides easy adjustment.
- Ultra low profile design eliminates height difference and facilitates natural gait.
- Fixed stirrups immobilise the ankle securely and comfortably.
- Impact resistant and shock absorbing sole for enhanced protection.
- Rocker Bottom sole allows easier ambulation.
- Padded soft foam keeps lower leg comfortable.
- Excellent durability

Fixed Walker Sizing Information:

Code	Size	UK Size	EU Size
FXPROS	Small	1-3.5	33-36
FXPROM	Medium	4-7	37-40
FXPROL	Large	7.5-9.5	41-44
FXPROXL	X-Large	10-14	45-48

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

INDICATIONS FOR USE:

- Replaces the need for heavy casts.
- Stable fracture of fibula, calcaneus and malleolus.
- Post ligament, soft tissue and tendon surgery.
- Post-surgery management.



Hook And Loop Closure



Fixed Stirrups



Rocker Bottom



Easier Ambulation



**INNER LINERS
ARE NOW
AVAILABLE**

BENECARE™ ROM WALKER

The **Benecare™ ROM Walker** is a high quality walker which conforms to the leg and immobilises the ankle with adjustable range of motion control.

FEATURES AND BENEFITS:

- Range of motion joint has 0-60 degrees in dorsiflexion and 0-70 degrees in plantar flexion adjustments in 10 degree increments.
- Hook & Loop closure provides easy adjustment.
- Ultra low profile design eliminates height difference and facilitates natural gait.
- Padded soft foam keeps lower leg comfortable.
- Metal stirrups immobilise the ankle securely and comfortably.
- Impact resistant and shock absorbing sole for enhanced protection.
- Rocker Bottom sole allows easier ambulation.

ROM Walker Sizing Information:

Code	Size	UK Size	EU Size
FXROMS	Small	1-3.5	33-36
FXROMM	Medium	4-7	37-40
FXROML	Large	7.5-9.5	41-44
FXROMXL	X-Large	10-14	45-48

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

INDICATIONS FOR USE:

- Replaces the need for serial casting.
- Stable fracture of fibula, calcaneus and malleolus.
- Post ligament, soft tissue and tendon surgery.
- Post-surgery Rehabilitation.
- The hinge allows controlled adjustments to the position of the ankle.
- Rupture of achilles tendon



Rocker Bottom



Flexion Adjuster



NEW

BENECARE™ ECO WALKER

The **BeneFoot™ Eco Walker** boot is an ultra-lightweight walker boot providing support and comfort to patients requiring short term immobilisation. The Eco Walker Boot comprises a rocker sole and removable foam lining helping improve short term rehabilitation capabilities



INDICATIONS FOR USE:

- Replaces the need for a heavy cast
- Stable fracture of fibula, calcaneus and malleolus
- Post ligament, soft tissue and tendon surgery
- Post-surgery management

FEATURES AND BENEFITS:

- Hook & Loop closure provides easy adjustment
- Ultra low profile design eliminates height difference and facilitates natural gait
- Fixed stirrups immobilise the ankle securely and comfortably
- Rocker bottom sole allows easier ambulation
- Padded soft foam keeps lower leg comfortable
- Ultra lightweight design makes it suitable for short term immobilisation

ROM Walker Sizing Information:			
Code	Size	UK Size	EU Size
FWS	Small	1-3,5	33-36
FWM	Medium	4-7	37-40
FWL	Large	7,5-9,5	41-44
FWXL	X-Large	10-14	45-48

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

41

BENECARE™ MEMORY FOAM ANKLE SUPPORT

The **Benecare™ Memory Foam** padding allows a custom fit ankle brace with lightweight foam providing a cushioned environment. The result is low profile custom fit comfort and custom length fitting. Also incorporated are two soft adjustable straps with Y hook clips, allowing custom length fitting.

INDICATIONS FOR USE:

- Post-operative stabilisation of the ankle
- Collateral ligament strains
- Inversion and eversion injuries
- Acute sport injuries
- Mild to moderate sprains
- Strains
- Post cast support
- Injury prevention

FEATURES AND BENEFITS:

- Two soft adjustable straps with loop.
- Removable Y hook clips allowing custom length fitting.
- Lightweight and durable.
- Quick and easy application for the clinician and patient.
- Adjustable distal heel strap.



Memory Foam



Adjustment Strap

Memory Ankle Brace Sizing Information:

Code	Size
MAB	Universal

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE



BENECARE™ AIR/GEL ANKLE BRACE

The BeneCare™ Air/Gel Ankle Brace has a removable Air top cell and a Gel lower cell that allows the brace to be used for cold therapy. Also incorporated are two soft adjustable straps with Y hook clips, allowing custom length fitting.

FEATURES AND BENEFITS:

- Two soft adjustable straps with loop.
- Removable Y hook clips allowing custom length fitting.
- Lightweight and breathable.
- Quicker application.



Gel Cushion Inner



Adjustable Straps

Air/Gel Ankle Brace Sizing Information:

Code	Size
AGAB	Universal Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

INDICATIONS FOR USE:

- Cold Therapy treatment
- Post-operative stabilisation of the ankle
- Collateral ligament strains
- Inversion and eversion injuries
- Acute sport injuries
- Mild to moderate sprains
- Strains
- Post cast support
- Injury prevention



BENECARE™ EXO ANKLE BRACE

The BeneCare™ EXO Ankle Brace is a superb lightweight ankle brace designed to give maximum protection to your ankles. The Exo Ankle Brace technology forms a plastic exoskeleton directly onto soft, breathable fabric, combining the compression and comfort of a soft ankle brace with the stability of a stirrup ankle brace.

FEATURES AND BENEFITS:

- The plastic exoskeleton frame minimizes inversion and eversion.
- Contours to fit ankle anatomy.
- Lace up strapping system for desired support.
- Lightweight fabric gives comfort and compliance.
- Ventilation holes allow the support to breathe and remain cool.
- Machine washable.

FITTING MEASUREMENTS:

- Ankle circumference diagonally over the base of the heel and across the flexor surface.

EXO Ankle Brace Sizing Information:

Code	Size	Ankle Circ. (cm)
EBABS	Small	305 - 330
EBABM	Medium	330 - 355
EBABL	Large	355 - 380
EBABXL	X-Large	380 - 405

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

INDICATIONS FOR USE:

- Ankle sprain
- Arthritis
- Contusion
- Tendonitis



BENECARE™ SILICON/GEL HEEL CUPS**FEATURES AND BENEFITS:**

- Uses dual durometer for maximum comfort while standing, walking or running.
- Medium-durometer silicone in main body efficiently absorbs impact and recovers quickly.
- Softer silicone blue zone under heel centre.
- Maximises shock absorption and comfort at point of impact.
- Retains shape - will not deform or degrade with use.
- Washable and reusable.

INDICATIONS FOR USE:

- Plantar fasciitis
- Heel spurs
- Heel pain
- Shock absorbent
- Bursitis
- Achilles tendonitis
- Pressure ulcers

**Silicone/Gel Heel Cups Sizing Information:**

Code	Size	UK Size	EU Size
SHC-S	Small	3-5	35-38
SHC-M	Medium	6-8	38.5-42.5
SHC-L	Large	9-11	43-45.5
SHC-XL	X-Large	12-14	46-49

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

**INDICATIONS FOR USE:**

- Post-op toe alignment.
- Hallux valgus.
- Hammer toe.
- Tailor's Bunion.

DARCO TOE ALIGNMENT SPLINT

The toe alignment splint is an effective way of obtaining excellent post operative toe alignment during the rehabilitation process.

FEATURES AND BENEFITS:

- New, thin elastic band fits easily into footwear and reduces slippage
- Toe straps are softer and will not irritate
- Special T-Strap that wraps comfortably around the great toe and fastens easily with hook and loop closures
- Soft foam loop straps slide easily over the nail to the base of the toe and is fixed to the metatarsal band
- Loop straps can be used after surgery to correct Tailor's bunions

Toe Alignment Splint Sizing Information:

Code	Size	Colour
TAS	Universal	White

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

ACHILLES HEEL PADS

The Benecare™ Achilles Heel Pad is designed to be inserted in the Airshell walking braces, the Heel Wedges provide additional Achilles tendon protection. The heel supports can be placed on top of each other to achieve the desired heel raise. Three heel supports equal 22° of heel raise. Two heel supports equal 16° of heel raise. One heel support equals 10° of heel raise

FEATURES AND BENEFITS:

- Adjustable height
- Achieve desired heel raise.
- Minimal patient discomfort

Achilles Heel Pads Sizing Information:

Code Left	Code Right	Pack
FXSWL	FXSWR	Pads x 5

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

