

ARM & SHOULDER

EXTENDER ARM BRACE.....	51
PRE-CUT MOULDED SPLINT.....	52
ELBOW CORK SPLINT	52
HUMERUS CORK SPLINT	52
SLING ON A ROLL.....	53
POLY ARM SLING.....	53





Adjustable Sliding Struts



Adjustable ROM Limitation



FEATURES AND BENEFITS

- Easy to use Extender hinge.
- Sliding struts lead to a more customised fit.
- Simple press and rotate ROM settings feature.
- Extension limitations: 0°, 15°, 30°, 45°, 60°, 75°, 90°.
- Flexion limitations: 0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°.
- Special left + right design.



BENECARE™ EXTENDER ARM BRACE

The **Benecare™ Extender Arm Brace** uses our proven technologies, the Extender elbow brace represents a new standard of design benefiting from telescopic struts. You will see improved patient outcome with superior comfort, clean aesthetics and simple hinge for post op and trauma indications.

INDICATIONS FOR USE

- Controlled elbow flexion/extension following surgery or trauma.
- Treatment of elbow dislocations.

Extender Arm Brace Sizing Information:

Code Left	Code Right	Size
EABL	EABR	Universal

PLEASE SEE PRICE LIST FOR NHS CODE

BENECAST™ PRE CUT MOULDED SPLINTS

The **BeneCast™ Pre-Cut Moulded Splints** are an excellent way of making a rigid splint quickly and efficiently with minimal mess. Simply wet the splint, by spraying or running under the tap and roll in a towel to absorb excess moisture. Apply directly to the limb using a Benecare cotton crepe bandage. The splint will mould perfectly to the contours of the limb and will reach initial set in 4 - 7 minutes and functional strength in 30 mins.

FEATURES AND BENEFITS:

- Can be used quickly and efficiently in A&E, theatres, plaster room, hand therapy, within the hospital setting and also by the ambulance service in emergency situations.
- 100% penetration of X-ray.
- 5-times lighter and stronger than Plaster Casting.
- Fast and efficient method to make a removable cast.
- Can be secured by either bandage or hook & loop tape.
- Pre-cut to its individual shape the splint requires minimal trimming or preparation.
- Excellent conformability to mould to contours of the limb.



Pre-Cut Moulded Arm Splints Sizing Information:

Code	Description
PCMS005	Long Arm

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE



Splint and hook & loop straps



Hook & loop straps

BENECARE™ ELBOW AND HUMERUS SPLINTS

The **UFS System** (Universal Fixation Splint) is an innovative, universal, ready to use splinting system made out of cork which is also the better substitute for Plaster of Paris (P.O.P.). As the system is an easy to use, fast and clean product with a high level of patient comfort it is a very suitable immobilisation product for both the A+E Department and the Fracture Clinic.

INSTRUCTIONS FOR USE

- Select the appropriate size of splint. If needed, the splint can easily be adjusted with scissors to the desired measurements
- Apply stockinette and under-cast padding suitable for the indication to protect the limb
- Secure to the limb with adhesive hook and loop, or Benecast cotton crepe bandage for additional compression.
- The UFS-splint will automatically follow the contours of the extremities



Elbow Splint Sizing Information:

Code	Description	Size
UFCS-04S	Elbow Splint Adult	Small
UFCS-04L	Elbow Splint Adult	Medium-Large

Paediatric Elbow Splint Sizing Information:

Code Left	Code Right	Size
UFCS-14S	Elbow Splint Paediatric	Small
UFCS-14M	Elbow Splint Paediatric	Medium
UFCS-14L	Elbow Splint Paediatric	Large

Humerus Splint Sizing Information:

Code Left	Code Right	Size
UFCS-06LS	UFCS-06RS	Small-Medium
UFCS-06LL	UFCS-06RL	Medium-Large

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

BENESLING™

The **BeneSling™** is a convenient immobilisation system for the upper extremity. The combination of a strap with the Y-tab hook closures offers users and patients an easy-to-adjust arm support.

FEATURES AND BENEFITS

- 24 metres of strapping
- Stretch resistant foam keeps arm in desired position
- Y-tab hook fasteners allow easy positioning adjustments of the arm
- Strap width remains constant during use which guarantees less compression in the neck
- 30 Y hook & loop tabs
- Custom fit for paediatric or adult with a single product

INDICATIONS FOR USE:

- Immobilisation of the upper extremity.



Y-tab hook fasteners



24 metres in box dispenser

BeneSling™ Sizing Information:

Code	Description	Size
BSR	2 rolls Universal fits left + right	5.5cm x 12m

PLEASE SEE PRICE LIST FOR NHS CODE



BENECARE™ POLY ARM SLING

The **BeneCare™ Poly Arm Sling** is a simple shoulder immobiliser to be used after surgery or trauma with the option to maintain the arm in raised position.

FEATURES AND BENEFITS:

- Maintains upper arm in raised position
- Breathable nylon/foam laminate construction
- Easy to apply with touch-and-close fastening
- Universal, left or right
- Washable at 40°C

INDICATIONS FOR USE:

- Post upper limb trauma
- Post upper limb surgery
- Hand surgery
- Fracture of the clavicle



Anti Twist Loop



Adjustable Waist Strap

BeneCare Poly Arm Sling Sizing Information:

Code	Description	Size
PAS	Arm Sling Universal fits left + right	Universal

PLEASE SEE PRICE LIST FOR NHS CODE

