


FIBER OPTIC  
**GREEN SYSTEM**  
DISPOSABLE BLADES




- Single patient use only.
- Latex free and safe in use.
- Especially very effective under emergency situations.
- Green system no direct contact of light.
- Stainless steel spatula and aluminum heel makes, practically unbreakable.
- Semi-metal blade also hold reinforcement metal hook to make it unbreakable.
- Low cost price with high quality.
- Round tip helps reduce trauma.
- Plastic optically fiber 5mm dia gives clear and maximum illumination.
- Black shrink tube on individual light guide protects light dispersion or reflection to doctor as well.
- Available in all sizes of Macintosh and Miller.
- Supplied surgically cleaned.
- Hygienically safe in use.
- Fully conforms to ISO 7376 and existing handles.
- Individually packed and labeled.
- In case of double polythene packing, reusable battery handle can be safe for use.
- Low profile shape of blades protects the patient from upper part of teeth.





Code	Description	Blade Number	Size
AL-1801-1	Infant	1	95mm
AL-1801-2	Child	2	105mm
AL-1801-3	Medium Adult	3	135mm
AL-1801-3.5	Ext. Medium Adult	3.5	145mm
AL-1801-4	Large Adult	4	155mm
AL-1801-5	X-Large	5	175mm



Code	Description	Blade Number	Size
AL-1802-00	Premature	00	65mm
AL-1802-0	Neonate	0	78mm
AL-1802-1	Infant	1	102mm
AL-1802-2	Child	2	154mm
AL-1802-3	Medium Adult	3	195mm
AL-1802-4	Large Adult	4	205mm



Code	Description	Blade Number	Size
AL-1803-1	Infant	1	95mm
AL-1803-2	Child	2	105mm
AL-1803-3	Medium Adult	3	135mm
AL-1803-3.5	Ext. Medium Adult	3.5	145mm
AL-1803-4	Large Adult	4	155mm
AL-1803-5	X-Large	5	175mm



Code	Description	Blade Number	Size
AL-1804-00	Premature	00	65mm
AL-1804-0	Neonate	0	78mm
AL-1804-1	Infant	1	102mm
AL-1804-2	Child	2	154mm
AL-1804-3	Medium Adult	3	195mm
AL-1804-4	Large Adult	4	205mm

## Fiber Optic Green System Disposable Blades

- Molded from highly rigid plastic material which makes it unbreakable
- Single patient use only
- Especially very effective under emergency situations.
- Green system no direct contact of light.
- Low cost price with high quality
- Round tip helps reduce trauma
- Plastic optically fiber 5mm dia gives clear and maximum illumination.
- Available in all sizes of Macintosh and Miller
- Hygienically safe in use
- Fully conforms to ISO 7376 and existing handles
- Individually packed and labeled

## Macintosh Blades, Fiber Optic Acrylic, Full-Plastic



Code	Description	Blade Number	Size
AL-1805-1	Infant	1	95mm
AL-1805-2	Child	2	115mm
AL-1805-3	Medium Adult	3	135mm
AL-1805-4	Large Adult	4	155mm

## Miller Blades Fiber Optic Acrylic, Full-Plastic



Code	Description	Blade Number	Size
AL-1806-0	Neonate	0	80mm
AL-1806-1	Infant	1	105mm
AL-1806-2	Child	2	155mm
AL-1806-3	Medium Adult	3	195mm
AL-1806-4	Large Adult	4	205mm

## MacFlex Blades, Fiber Optic Acrylic, Full-Metalic



- Very useful while a difficult intubation is encountered.
- Less force required to view the larynx as compare to normal Macintosh blades.
- Single patient use only.
- Especially very effective during minimal neck movement of the patient.
- Stainless steel spatula and aluminum heel makes, practically unbreakable.
- Low cost price with high quality.
- Round tip helps reduce trauma.
- Plastic optically fiber 5mm dia gives clear and maximum illumination.
- Supplied surgically cleaned
- Hygienically safe in use
- Fully conforms to ISO 7376 and existing handles
- Individually packed and labeled
- Low profile shape of blades protects the patient from upper part of teeth.

Code	Description	Blade Number	Size
AL-1807-2	Child	2	105mm
AL-1807-3	Medium Adult	3	135mm
AL-1807-4	Large Adult	4	155mm

- Knurled barrel gives strong grip to the user.
- Single patient use only
- LED light source with 25000 LUX
- LED handles with built in batteries
- Available in fullmetal and semimetal construction.
- Very economical in pricing with high quality
- Crafted from light weight aluminium metal.
- Supplied surgically clean which makes it hygienically safe in use.
- Fully conforms to ISO 7376 and existing blades as well.
- Individually packed and labeled.

### Battery Handles Fiber Optic LED, Semimetal with Batteries

Code	Description	Battery	Size
AL-1851-2	Infant		115mm
AL-1851-3	Pediatric		150mm
AL-1851-4	Medium		150mm

### Battery Handles Fiber Optic LED, Metallic with Batteries

Code	Description	Battery	Size
AL-1852-2	Infant		115mm
AL-1852-3	Pediatric		150mm
AL-1852-4	Medium		150mm

### Battery Handles Fiber Optic LED, Semimetal without Batteries

Code	Description	Battery	Size
AL-1858-2	Infant		115mm
AL-1858-3	Pediatric		150mm
AL-1858-4	Medium		150mm

### Battery Handles Fiber Optic LED, Metallic without Batteries

Code	Description	Battery	Size
AL-1859-2	Infant		115mm
AL-1859-3	Pediatric		150mm
AL-1859-4	Medium		150mm

Disposable Metallic Kits (Battery Included)



Code	Description
AL-1853-00	Miller Blade 00 with Battery Handle Med
AL-1853-0	Miller Blade 0 with Battery Handle Med
AL-1853-1	Miller Blade 1 with Battery Handle Med
AL-1853-2	Miller Blade 2 with Battery Handle Med
AL-1853-3	Miller Blade 3 with Battery Handle Med
AL-1853-4	Miller Blade 4 with Battery Handle Med

Code	Description
AL-1854-0	Mac Blade 0 with Battery Handle Med
AL-1854-1	Mac Blade 1 with Battery Handle Med
AL-1854-2	Mac Blade 2 with Battery Handle Med
AL-1854-3	Mac Blade 3 with Battery Handle Med
AL-1854-4	Mac Blade 4 with Battery Handle Med
AL-1854-5	Mac Blade 5 with Battery Handle Med

Disposable Metallic Kits (W/O Batteries)



Code	Description
AL-1855-00	Miller Blade 00 with Battery Handle Med
AL-1855-0	Miller Blade 0 with Battery Handle Med
AL-1855-1	Miller Blade 1 with Battery Handle Med
AL-1855-2	Miller Blade 2 with Battery Handle Med
AL-1855-3	Miller Blade 3 with Battery Handle Med
AL-1855-4	Miller Blade 4 with Battery Handle Med

Code	Description
AL-1856-0	Mac Blade 0+Battery Handle Med
AL-1856-1	Mac Blade 1+Battery Handle Med
AL-1856-2	Mac Blade 2+Battery Handle Med
AL-1856-3	Mac Blade 3+Battery Handle Med
AL-1856-4	Mac Blade 4+Battery Handle Med
AL-1856-5	Mac Blade 5+Battery Handle Med

Replacement LED Capsule

Code	Description
AL-1857-1	Replacement LED capsule