


CONVENTIONAL
DISPOSABLE
BLADES




- Single patient use only.
- Cool and bright illumination as compare to conventional vacuum lamp.
- Latex free and safe in use.
- Especially very effective under emergency situations.
- LED Light source with 5000 LUX at the tip of the blade.
- Stainless steel spatula and aluminum heel makes, practically unbreakable.
- Low cost price with high quality.
- Round tip helps reduce trauma.
- Fully conforms to ISO 7376 Standard.
- Available in all sizes of Macintosh and Miller.
- Supplied surgically cleaned.




Code	Description	Blade Number	Size
AL-1601-1	Infant	1	95mm
AL-1601-2	Child	2	105mm
AL-1601-3	Medium Adult	3	135mm
AL-1601-4	Large Adult	4	155mm
AL-1601-5	X-Large	5	175mm



Code	Description	Blade Number	Size
AL-1602-00	Premature	00	65mm
AL-1602-0	Neonate	0	78mm
AL-1602-1	Infant	1	102mm
AL-1602-2	Child	2	154mm
AL-1602-3	Medium Adult	3	195mm
AL-1602-4	Large Adult	4	205mm



Code	Description	Blade Number	Size
AL-1603-1	Infant	1	95mm
AL-1603-2	Child	2	105mm
AL-1603-3	Medium Adult	3	135mm
AL-1603-3.5	Ext. Medium Adult	3.5	145mm
AL-1603-4	Large Adult	4	155mm
AL-1603-5	X-Large	5	175mm



Code	Description	Blade Number	Size
AL-1604-00	Premature	00	65mm
AL-1604-0	Neonate	0	78mm
AL-1604-1	Infant	1	102mm
AL-1604-2	Child	2	154mm
AL-1604-3	Medium Adult	3	195mm
AL-1604-4	Large Adult	4	205mm

Battery Handles for Conventional Disposable Blades

- Knurled barrel gives strong grip to the user.
- Single patient use only.
- Compatible LED blades only.
- No need to install batteries because batteries are included.
- Available in full-metal and semimetal construction.
- Very economical in pricing with high quality
- Crafted from light weight aluminium metal.
- Supplied surgically clean which makes it hygienically safe in use.
- Fully conforms to ISO 7376 and existing blades as well.



Battery Handles Conventional, Semimetal

Code	Description	Battery	Size
AL-1651-2	Infant		115mm
AL-1651-3	Pediatric		150mm
AL-1651-4	Medium		150mm



Battery Handles Conventional, Metallic

Code	Description	Battery	Size
AL-1652-1	Infant,		115mm
AL-1652-3	Pediatric		150mm
AL-1652-4	Medium		150mm

