



KNEE

HINGED KNEE STABILISER.....	45
NEO-WRAP KNEE BRACE.....	46
LONG NEO-WRAP BRACE.....	46
ROM KNEE BRACE.....	47
TRI PANEL KNEE SPLINT.....	48
PRO-ACTION KNEE BRACE	49

BENECARE™ HINGED KNEE STABILISER

This **Benecare™ Hinge Knee Stabiliser** offers support for weak, unstable knees, and protection during rehabilitation.

FEATURES AND BENEFITS:

- Polycentric hinge helps prevent hyperextension of the knee without restricting mobility.
- Soft neoprene blend material provides soothing warmth by retaining body heat, keeping the knee joint and soft tissues flexible.
- Opposing adjustable front and back straps hold brace in place.
- Contoured construction for comfortable fit minimizes bulk behind knee joint.
- Adjustable straps provide controlled compression for custom fit.
- Patella buttress stabilises and protects kneecap.
- Increasing medial / lateral support.

INDICATIONS FOR USE:

- Rheumatoid arthritis
- Mild to moderate hyperextension
- Mild to moderate valgus/ varus deformity

FITTING INFORMATION:

- Circumference at mid-patella.

Hinged Knee Stabiliser Sizing Information:

Code	Size	Knee circ. (cm)
HKSS	Small	33-37
HKSM	Medium	37-41
HKSL	Large	41-45
HKSXL	X-Large	45-50
HKSSXL	XX-Large	51-55

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE





Short Neo-Wrap Brace



Long Neo-Wrap Brace

BENECARE™ NEO-WRAP BRACE

The **Benecare™ Neo-Wrap Brace** front opening design and easy-opening tabs on the Velcro™ straps make this brace suitable for the arthritic or elderly patient, who may suffer from limited dexterity. The brace incorporates hyper stop hinges provide medial and lateral support for unstable knee, while allowing flexion. The internal access to the hinge allows for easy removal of the hinge for laundering.

FEATURES AND BENEFITS:

- Front opening.
- Universal left & right.
- Contrasting proximal & distal strap to counter rotation.
- Open popliteal area for added comfort during flexion.
- Removable hyper stop hinges.
- Easy grip straps for patients with poor dexterity.
- Fully adjustable front & rear straps

INDICATIONS FOR USE:

- Increasing medial/lateral support
- Rheumatoid arthritis
- Mild to moderate hyperextension
- Mild to moderate valgus/varus deformity



Removable Hinges

Neo-Wrap Brace Sizing Information:

Code-Beige	Size	Knee circ. (cm)	Code-Black	Size	Knee circ. (cm)
SKKBXS	X-Small	25 - 28	BSKKBXS	X-Small	25 - 28
SKKBS	Small	28 - 34	BSKKBBS	Small	28 - 34
SKKBM	Medium	34 - 40	BSKKBM	Medium	34 - 40
SKKBL	Large	40 - 46	BSKKBBL	Large	40 - 46
SKKBXL	X-Large	46 - 56	BSKKBXL	X-Large	46 - 56
SKKBXXL	XX-Large	56 - 61	BSKKBXXL	XX-Large	56 - 61

Long Neo-Wrap Brace Sizing Information:

Code-Beige	Size	Knee circ. (cm)	Code-Black	Size	Knee circ. (cm)
LSKKB-XS	X-Small	25 - 28	BLSKKB-XS	X-Small	25 - 28
LSKKB-S	Small	28 - 34	BLSKKB-S	Small	28 - 34
LSKKB-M	Medium	34 - 40	BLSKKB-M	Medium	34 - 40
LSKKB-L	Large	40 - 46	BLSKKB-L	Large	40 - 46
LSKKB-XL	X-Large	46 - 56	BLSKKB-XL	X-Large	46 - 56
LSKKB-XXL	XX-Large	56 - 61	BLSKKB-XXL	XX-Large	56 - 61

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

BENECARE™ ROM KNEE BRACE

The **BeneCare™ ROM Knee Brace** provides excellent collateral stabilisation via long side-bars and prevents valgus or varus deviation as well as preventing hyperextension using the braces adjustable dial lock hinges.

FEATURES AND BENEFITS:

- Dial lock hinges with adjustable ROM limitation
- Extension limitation at: 0°, 15°, 30°, 45°, 60°, 75°
- Flexion limitation at: 0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
- Immobilisation limitation at: 0°, 15°, 30°, 45°, 60°, 75°
- Laminated soft foam material provides excellent conformability to the contours of the limb offering additional patient comfort
- Length: 18", 20", 22", 24"

INDICATIONS FOR USE:

- Stable and unstable knee injuries
- Post- operative immobilisation for
 - ACL/PCL
 - Collateral Ligaments
 - Tibial Plateau Fractures

ROM Knee Brace Sizing Information:		
Code	Size	Length
RKB16	Small	16" (41cm)
RKB20	Medium	20" (51cm)
RKB24	Large	24" (61cm)

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE



Adjustable ROM Limitation



Adjustable Straps



BENECARE™ TRI PANEL KNEE SPLINT



Tri Panel Knee Splint Sizing Information:

Code	Size	Length (inches)
TPKS10	XX-Small	10" (25.5cm)
TPKS12	X-Small	12" (30.5cm)
TPKS14	Small	14" (36cm)
TPKS16	Medium	16" (41cm)
TPKS18	Large	18" (46cm)
TPKS20	X-Large	20" (51cm)
TPKS22	XX-Large	22" (56cm)
TPKS24	XXX-Large	24" (61cm)

Universal fits left + right

PLEASE SEE PRICE LIST FOR NHS CODE

BENECARE™ TRI PANEL KNEE SPLINT

The BeneCare™ Tri-Panel Knee Splint offers patients a high quality support designed for short and long term use. The detachable panels can be adjusted to ensure a perfect fit.

FEATURES AND BENEFITS:

- Three piece design using high quality foam for a precise fit.
- Removable medial, lateral and posterior stays.
- Comfortable, washable, skin-friendly terry towelling cloth material.
- Hook & loop straps with loop closures for simple fitting.

INDICATIONS FOR USE:

- Post-operative immobilisation.
- Pre-operative immobilisation.
- Patella dislocation.
- Conservative treatment of ligament injury.
- Acute anterior knee pain.

BENECARE™ PRO ACTION KNEE

The **BeneCare™ Pro Action Knee Brace** has a high quality polycentric hinge design that fits the actual movement of the human knee, facilitating enhanced rehabilitation and support. Its rigid construction prevents hyperextension via the adjustable ROM hinge and it contains a 4 way leverage design for optimal patient stabilisation.

FEATURES AND BENEFITS:

- ROM hinges with plug-in flex/ext stops.
- Extension limitation at: 0°, 10°, 20°, 30°.
- Flexion limitation at: 0°, 30°, 45°, 60°, 90°.
- Immobilisation limitation at: 0°, 30°.
- Hypoallergenic silicone-grip strap pads for anti-migration.
- Concentric condyle pads correspond to bilateral femoral condyles.
- Numbered straps provide desirable counter force and adjust optimal size.

INDICATIONS FOR USE:

- Knee instability due to anterior and/or posterior cruciate ligament (ACL/PCL) injury/rupture.
- To restrict the range of motion after reconstructive surgery.
- Recommended for low to medium impact sports.

FITTING INFORMATION:

- Thigh Sizing Guide: measure 15cm (6") above the medial condyle.
- Calf Sizing Guide: measure 12cm (5") below the medial condyle.



Hinges



Hinges



Hypoallergenic Straps



Immobilisation Limitation

Pro Action Knee Brace Sizing Information:

Code Left	Code Right	Size	Thigh Circ.	Calf circ
PAKBLS	PAKBRS	Small	15.5" - 18.5" (39-47cm)	11.5" - 14.5" (29-37cm)
PAKBLM	PAKBRM	Medium	18.5" - 21" (47-53.5cm)	14.5" - 17" (37-43cm)
PAKBLL	PAKBLL	Large	21" - 23.5" (53.5-60cm)	17" - 19.3" (43-50cm)
PAKBLXL	PAKBRXL	X-Large	23.5" - 26.5" (60-67cm)	19.3" - 22.5" (50-57cm)

PLEASE SEE PRICE LIST FOR NHS CODE