

Video Laryngoscope UED-C series



UED-C



- Infection control
- Cost control

ZHEJIANG UE MEDICAL CORP.

✉ EXPORT@UEWORLD.COM

☎ +86 571 28057363

🌐 WWW.UEWORLD.COM/EN.PHP

📍 8 YOUYI RD, BAITA ECONOMIC
DEVELOPMENT ZONE, XIANJU,
ZHEJIANG 317300, CHINA

UED-M3-S adds on:



TOUCH SCREEN

- image mirroring: flips upside-down whenever needed
- auto shutdown: powers off after 10 mins of inactivity

WIFI READY

- real-time screen sharing
- fast access to photos and videos

2 monitor options (UED-M3, UED-M3-S)



UED-C4
adult, obesity



UED-C3
6-adult



UED-C2
2-6 yrs



UED-MC1
0-2 yrs



UED-C1
0-2 yrs



UED-C0
newborn



Compatible with UE M8 8" monitor

MAINTAINED QUALITY

- reusable slim handle
- signature blades design
- anti-fog coating

CLEARER VIEW

- upgraded monitor
- advanced camera

MORE FOR CLINICAL NEEDS

- for all ages and sizes
- a compatible monitor for a series of products

SPECIFICATIONS

MONITOR		BLADES	
Display	3" TFT LCD, 720x480 px	Angle of View	60°
Tilt and Rotation	0°-110° tilt, 0°-270° rotate	Light Source	high intensity LED
Memory	32G, take > 400k photos, record ≥16 hrs of video		
Battery Life	> 3.5 hrs on a full charge		
Charging Time	< 3 hrs		