

# OUR ACCESSORIES

Various options are available to make the system adaptable to any situation.



Nocospray is fully customizable.



## Remote control

The remote control plug enables the launch of Nocospray from outside the room.



## Timer

Plugged into both a mains socket and Nocospray, the timer enables the delay or repeat treatments.



## Spraying gun

Attached to a Nocospray thanks to a flexible hose, the spraying gun treats areas difficult to access (ventilation ducts...)



## Wall-crossing system

To disinfect a room without entering it, the wall-crossing system is particularly adapted to vehicles and white rooms.



## Wall mounted Nocospray IS

Same features as Nocospray but fixed on the wall of the room to be treated.



## Wall mounted Nocospray OS

Fixed on the wall, it enables the treatment of a room from the outside thanks to a wall-crossing system.



## Double nozzle

The double nozzle enables a bidirectional diffusion specially adapted to long and narrow places such as corridors.



## Traceability

The traceability system connects the Nocospray to a computer thanks to a USB cable, to record all treatments.

## CONTACT



**DIRECT STEP**  
الخطوة المباشرة

Direct Step Co. for Specialized Service

Tel: +966 12 652 55 49

Fax: +966 12 652 55 48

Kingdom of Saudi Arabia

P.O. Box 7769, Jeddah 21472