

SaCoVLM™

Seeing is believing

Disposable Video Laryngeal Mask

Product Specification

Model	Specification	Patient Weight (kg)	Largest ETT (mm)	Suitable drainage tube (Fr)	Instrument for Airway Tube, Max. OG Size (mm)
MLM30C	3.0#	30-50	7.0	16	5.8
MLM40C	4.0#	50-70	7.5	16	5.8
MLM50C	5.0#	70-100	8.0	16	5.8
DLM30	3.0#	30-50	7.0	16	5.8
DLM40	4.0#	50-70	7.5	16	5.8
DLM50	5.0#	70-100	8.0	16	5.8



SaCoVLM™

Seeing is believing

Disposable Video Laryngeal Mask

ZHEJIANG UE MEDICAL CORP.

✉ EXPORT@UEWORLD.COM

🌐 WWW.UEWORLD.COM/EN.PHP

☎ +86 571 28057363

📍 8 YOUYI RD, BAITA ECONOMIC DEVELOPMENT ZONE, XIANJU, ZHEJIANG 317300, CHINA

SaCoVLM™

Seeing is believing

Disposable Video Laryngeal Mask

Clinical Uses

1. Optimal for emergent airways
2. Rapid ventilation for normal and difficult airways
3. The preferred choice for rescue ventilation in various clinical scenarios
4. Video-assisted tracheal intubation and extubation
5. Video-assisted clinical diagnosis and treatment

Cleaning line

Inflation line

Visual tube channel

Ventilation (intubation) channel

Gastric tube channel

Pre-molded


ergonomic & for easier insertion and placement

Unique sloping distal end

facilitates accurate insertion of ETT

Inflatable cuff

large contact area lowers mucosal pressure & higher airtight seal avoids aspiration



- **Seeing is believing**
The visible insertion process allows accurate placement and therefore avoids air leakage.
- **Intraoperative monitoring**
Monitors the glottis for any case of reflux to reduce aspiration.
- **Facilitates safe intubation**
Guides the tracheal intubation under visible conditions, with safety and comfort.