

Video Laryngoscope with Single-use Blade



UED-D



- **Maximize efficiency**
- **Minimize cross-contamination**

ZHEJIANG UE MEDICAL CORP.

✉ EXPORT@UEWORLD.COM

☎ +86 571 28057363

🌐 WWW.UEWORLD.COM/EN.PHP

📍 8 YOUYI RD, BAITA ECONOMIC
DEVELOPMENT ZONE, XIANJU,
ZHEJIANG 317300, CHINA



UED-M3-S adds on:



TOUCH SCREEN

- image mirroring: flips upside-down whenever needed
- auto shutdown: powers off after 10 mins of inactivity



WIFI READY

- real-time screen sharing
- fast access to photos and videos

2 monitor options (UED-M3, UED-M3-S)



UED-D0

newborn, 2.7-3.8kg



UED-D1

0-2 yrs, 3.4-13.2kg



UED-D2

2-6 yrs, 13.2-23.6kg



UED-D3

6-adult, above 23.6kg



UED-D4

adult, obesity

SINGLE-USE BLADES

- designed to prohibit any reuse

GRAB AND GO

- quick set up
- go-to device for prehospital (EMS) & hospital
- no concerns about maintenance, cleaning, and disinfecting

WHENEVER YOU WANT IT

- lighter
- longer powered-up
- faster blades changing (hot-swappable)

SPECIFICATIONS

MONITOR		BLADES	
Display	3" TFT LCD, 720x480 px	Camera	CMOS
Tilt and Rotation	0°-110° tilt, 0°-270° rotate	Angle of View	60°
Memory	32G, take > 400k photos, record ≥16 hrs of video	Light Source	high intensity LED
Weight	150g	Weight	40g
Battery Life	> 3.5 hrs on a full charge		
Charging Time	< 3 hrs		