

# Reusable Video Laryngoscope



VL310



VL300



- **Cost Effective**
- **For patients of all ages and sizes**

ZHEJIANG UE MEDICAL CORP.

✉ EXPORT@UEWORLD.COM

☎ +86 571 28057363

🌐 WWW.UEWORLD.COM/EN.PHP

📍 8 YOUYI RD, BAITA ECONOMIC DEVELOPMENT ZONE, XIANJU, ZHEJIANG 317300, CHINA



## ALL-ANGLE ADJUSTABLE MONITOR

- VL310 additional key: take & store photos and videos easy to teach

## SMALL MONITOR ON A STREAMLINED SLIM HANDLE

- ergonomic & portable
- one hand operation
- a single person is sufficient for airway management

## SIGNATURE BLADES DESIGN

- visualize the glottis > fast & efficient
- hide the epiglottis > less vagal stimulation & safe

## OFFER EACH SIZE WITH CARE AND CONSIDERATION

- not simply scaling
- for all ages and sizes



**310-0 / 300SS**  
newborn



**310-1 / 300S**  
0-2yrs



**310-2 / 300SM**  
2-6yrs



**310-3 / 300M**  
6yrs+



**310-4 / 300L**  
adult, obesity

## SPECIFICATIONS

Display	2.7" TFT LCD	Battery Life	> 120 mins on a full charge
Tilt and Rotation	0°-130° tilt, 0°-270° rotate	Charging Time	< 3 hrs
Angle of View	60°	Memory	32G, take > 400k photos record ≥16 hrs of video